

## Peptide Array, Hepatitis C Virus, H77, NS4B Protein

### Catalog No. NR-3754

This reagent is the tangible property of the U.S. Government.

### For research use only. Not for human use.

#### Contributor:

BEI Resources

#### Manufacturer:

Bio-Synthesis, Inc.

#### Product Description:

The 40-peptide array spans the NS4B protein of hepatitis C virus, H77 (genotype 1a; GenPept: AAB67036).<sup>1</sup> Peptides are 12- to 18-mers, with 11 or 12 amino acid overlaps. Please see Table 1 for length and sequence of individual peptides.

#### Material Provided:

Peptides are provided lyophilized at 1 mg per vial.

#### Packaging/Storage:

Lyophilized peptides should be placed in a closed dry environment with desiccants and stored at -20°C or colder immediately upon arrival. A frost-free freezer should be avoided, since changes in moisture and temperature may affect peptide stability.

#### Solubility:

Solubility may vary based on the amino acid content of the individual peptide (see Table 2).

#### Reconstitution:

Lyophilized peptides should be warmed to room temperature for 1 hour prior to reconstitution. They should be dissolved at the highest possible concentration, and then diluted with water or buffer to the working concentration. Buffer should be added only after the peptide is completely in solution because salts may cause aggregation.

The most common dissolution process is 1 mg of peptide in 1 mL of sterile, distilled water. Peptides that are not soluble in water can almost always be dissolved in DMSO. Once a peptide is in solution, the DMSO can be slowly diluted with aqueous medium. Care must be taken to ensure that the peptide does not begin to precipitate out of solution. For cell-based assays, 0.5% DMSO in medium is usually well-tolerated.

Sonication and/or the addition of small amounts of dilute (10%) aqueous acetic acid for basic peptides, aqueous ammonia for acidic peptides or acetonitrile may also help dissolution (see Table 2). These solvents may not be

appropriate for certain applications, including cell-based assays.

#### Storage of Reconstituted Peptides:

The shelf life of peptides in solution is very limited, especially for sequences containing cysteine, methionine, tryptophan, asparagine, glutamine, and N-terminal glutamic acid. In general, peptides may be aliquoted and stored in solution for a few days at -20°C or colder. For long-term storage, peptides should be re-lyophilized and stored at -20°C or colder. If long-term storage in solution is unavoidable, peptide solutions should be buffered to pH 5-6, aliquoted and stored at -20°C or colder. Freeze-thaw cycles should be avoided.

#### Citation:

Acknowledgment for publications should read "The following reagent was obtained through BEI Resources, NIAID, NIH: Peptide Array, Hepatitis C Virus, H77, NS4B Protein, NR-3754."

#### Biosafety Level: 1

Appropriate safety procedures should always be used with this material. Laboratory safety is discussed in the following publication: U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention, and National Institutes of Health. Biosafety in Microbiological and Biomedical Laboratories. 5th ed. Washington, DC: U.S. Government Printing Office, 2007; see [www.cdc.gov/od/ohs/biosfty/bmb15/bmb15toc.htm](http://www.cdc.gov/od/ohs/biosfty/bmb15/bmb15toc.htm).

#### Disclaimers:

You are authorized to use this product for research use only. It is not intended for human use.

Use of this product is subject to the terms and conditions of the BEI Resources Material Transfer Agreement (MTA). The MTA is available on our Web site at [www.beiresources.org](http://www.beiresources.org).

While BEI Resources uses reasonable efforts to include accurate and up-to-date information on this product sheet, neither ATCC® nor the U.S. Government make any warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. Neither ATCC® nor the U.S. Government warrants that such information has been confirmed to be accurate.

This product is sent with the condition that you are responsible for its safe storage, handling, use and disposal. ATCC® and the U.S. Government are not liable for any damages or injuries arising from receipt and/or use of this product. While reasonable effort is made to ensure authenticity and reliability of materials on deposit, the U.S. Government, ATCC®, their suppliers and contributors to BEI Resources are not liable for damages arising from the misidentification or misrepresentation of products.

**Use Restrictions:**

**This material is distributed for internal research, non-commercial purposes only.** This material, its product or its derivatives may not be distributed to third parties. Except as performed under a U.S. Government contract, individuals contemplating commercial use of the material, its products or its derivatives must contact the contributor to determine if a license is required. U.S. Government contractors may need a license before first commercial sale. This material may be subject to third party patent rights.

**References:**

1. Yanagi, M., et al. "Transcripts from a Single Full-length cDNA Clone of Hepatitis C Virus Are Infectious When Directly Transfected into the Liver of a Chimpanzee." *Proc. Natl. Acad. Sci. U. S. A.* 94 (1997): 8738-8743. PubMed: 9238047. GenPept: AAB67036.

ATCC® is a trademark of the American Type Culture Collection.



Table 1		
Peptide	Length	Sequence
1 of 40	18	1 SQHLPYIEQGMMLAEQFK 18
2 of 40	18	8 EQGMMLAEQFKQKALGLL 25
3 of 40	18	15 EQFKQKALGLLQTASRHA 32
4 of 40	18	22 LGLLQTASRHAEVITPAV 39
5 of 40	18	29 SRHAEVITPAVQTNWQKL 46
6 of 40	18	36 TPAVQTNWQKLEVFVWAKH 53
7 of 40	16	43 WQKLEVFVWAKHMMWNFI 58
8 of 40	18	48 VFVWAKHMMWNFISGIQYLA 65
9 of 40	16	55 WNFISGIQYLAGLSTL 70
10 of 40	18	60 GIQYLAGLSTLPGNPAIA 77
11 of 40	18	67 LSTLPGNPAIASLMAFTA 84
12 of 40	17	74 PAIASLMAFTAAVTSPL 90
13 of 40	18	80 MAFTAAVTSPLTTGQTLL 97
14 of 40	18	87 TSPLTTGQTLLFNILGGW 104
15 of 40	18	93 GQTLLFNILGGWVAAQLA 110
16 of 40	18	100 ILGGWVAAQLAAPGAATA 117
17 of 40	18	107 AQLAAPGAATAFVVGAGLA 124
18 of 40	18	114 AATAFVVGAGLAGAAIGSV 131
19 of 40	18	121 AGLAGAAIGSVGLGKVLV 138
20 of 40	17	128 IGSVGLGKVLVDILAGY 144
21 of 40	18	134 GKVLVDILAGYGAGVAGA 151
22 of 40	18	141 LAGYGAGVAGALVAFKIM 158
23 of 40	15	148 VAGALVAFKIMSGEV 162
24 of 40	18	152 LVAFKIMSGEVPSTEDLV 169
25 of 40	18	159 SGEVPSTEDLVNLLPAIL 176
26 of 40	18	166 EDLVNLLPAILSPGALVV 183
27 of 40	18	173 PAILSPGALVVGVVCAAI 190
28 of 40	16	180 ALVVGVVCAAILRRHV 195

Table 1		
Peptide	Length	Sequence
29 of 40	18	185 VVCAAILRRHVGPGEAV 202
30 of 40	18	192 RRHVGPGEAVQWMNRLI 209
31 of 40	16	199 EGAVQWMNRLIAFASR 214
32 of 40	15	204 WMNRLIAFASRGNHV 218
33 of 40	17	208 LIAFASRGNHVSPHYV 224
34 of 40	18	214 RGNHVSPHYVPESDAAA 231
35 of 40	17	221 THYVPESDAAARVTAIL 237
36 of 40	16	227 SDAAAR/TAISSLTV 242
37 of 40	18	232 RVTAILSSLTVTQLLRR 249
38 of 40	15	239 SLTVTQLLRRRLHQWI 253
39 of 40	18	243 TQLLRRRLHQWISSECTP 260
40 of 40	12	250 HQWISSECTTPC 261

Table 2		
Peptide	Solubility	Solvent
1 of 40	1 mg/mL	6 M guanidine-HCl
2 of 40	1 mg/mL	100% DMSO
3 of 40	1 mg/mL	6 M guanidine-HCl
4 of 40	1 mg/mL	6 M guanidine-HCl
5 of 40	1 mg/mL	6 M guanidine-HCl
6 of 40	1 mg/mL	100% DMSO
7 of 40	1 mg/mL	100% DMSO
8 of 40	1 mg/mL	100% DMSO
9 of 40	1 mg/mL	100% DMSO
10 of 40	1 mg/mL	6 M guanidine-HCl
11 of 40	1 mg/mL	100% DMSO
12 of 40	1 mg/mL	100% DMSO
13 of 40	1 mg/mL	100% DMSO
14 of 40	1 mg/mL	100% DMSO
15 of 40	1 mg/mL	100% DMSO
16 of 40	1 mg/mL	100% DMSO
17 of 40	1 mg/mL	100% DMSO
18 of 40	1 mg/mL	100% DMSO
19 of 40	1 mg/mL	100% DMSO
20 of 40	1 mg/mL	100% DMSO
21 of 40	1 mg/mL	100% DMSO

Table 2		
Peptide	Solubility	Solvent
22 of 40	1 mg/mL	100% DMSO
23 of 40	1 mg/mL	100% DMSO
24 of 40	1 mg/mL	100% DMSO
25 of 40	1 mg/mL	6 M guanidine-HCl
26 of 40	1 mg/mL	100% DMSO
27 of 40	1 mg/mL	100% DMSO
28 of 40	1 mg/mL	100% DMSO
29 of 40	1 mg/mL	6 M guanidine-HCl
30 of 40	1 mg/mL	6 M guanidine-HCl
31 of 40	1 mg/mL	100% DMSO
32 of 40	1 mg/mL	100% DMSO
33 of 40	1 mg/mL	6 M guanidine-HCl
34 of 40	1 mg/mL	70% acetonitrile in water
35 of 40	1 mg/mL	6 M guanidine-HCl
36 of 40	1 mg/mL	100% DMSO
37 of 40	1 mg/mL	100% DMSO
38 of 40	1 mg/mL	100% DMSO
39 of 40	1 mg/mL	100% DMSO
40 of 40	1 mg/mL	70% acetonitrile in water