

**SARS-Related Coronavirus 2, Wuhan-Hu-1 Spike D614G-Pseudotyped Lentivirus, Luc2/ZsGreen**

**Catalog No. NR-53819**

**For research use only. Not for use in humans.**

**Contributor:**

Jesse Bloom, Associate Professor, Department of Basic Sciences, Fred Hutchinson Research Center, Seattle, Washington, USA and Alejandro Balasz, Assistant Professor, The Ragon Institute of MGH, MIT and Harvard, Cambridge, Massachusetts, USA

**Manufacturer:**

BEI Resources

**Product Description:**

A pseudotyped lentiviral form of the spike (S) glycoprotein gene from severe acute respiratory syndrome-related coronavirus 2 (SARS-CoV-2), Wuhan-Hu-1 (GenBank: [NC\\_045512](#)) was produced in human embryonic kidney HEK293 cells and purified by sucrose cushion. NR-53819 expresses a C-terminally truncated S glycoprotein containing a D614G mutation. The mutation increases titers of viral particles pseudotyped with SARS-CoV-2 S glycoprotein, as well as synthetic firefly luciferase (Luc2) and synthetic *Zoanthus* sp. green fluorescent protein (ZsGreen1).<sup>1,2</sup> NR-53819 is produced from a plasmid kit (BEI Resources NR-52516, NR-52517, NR-52518, NR-52519 and NR-53765; kit NR-53817). Protocols for the use of these items are published, and updates are available at [protocols.io](#).<sup>1,3</sup>

The S glycoprotein mediates viral binding to the host angiotensin converting enzyme 2 (ACE2). This protein forms a trimer, and when bound to a host receptor, allows fusion of the viral and cellular membranes. The S protein is a target for neutralizing antibodies.<sup>4</sup>

See Appendix I for assay information. Successful use of this product necessarily depends on third party reagents. It is critically important that users confirm effective assay performance in their own laboratories using the appropriate control templates.

**Material Provided:**

Each vial contains approximately 10 µL of purified lentiviral particles from transfected HEK293 supernatant in DMEM supplemented with 10% heat-inactivated fetal bovine serum. NR-53819 includes sufficient material for one 96 well plate neutralization assay.

**Packaging/Storage:**

NR-53819 was packaged aseptically in screw-capped plastic cryovials. The product is provided frozen and should be stored at -60°C or colder immediately upon arrival. Freeze-thaw cycles should be avoided.

**Citation:**

Acknowledgment for publications should read “The following reagent was obtained through BEI Resources, NIAID, NIH: SARS-Related Coronavirus 2, Wuhan-Hu-1 Spike D614G-Pseudotyped Lentivirus, Luc2/ZsGreen, NR-53819.”

**Biosafety Level: 2**

Appropriate safety procedures should always be used with this material. Laboratory safety is discussed in the following publication: U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention, and National Institutes of Health. [Biosafety in Microbiological and Biomedical Laboratories](#). 6th ed. Washington, DC: U.S. Government Printing Office, 2020; see [www.cdc.gov/biosafety/publications/bmbli5/index.htm](http://www.cdc.gov/biosafety/publications/bmbli5/index.htm).

**Disclaimers:**

You are authorized to use this product for research use only. It is not intended for human use.

Use of this product is subject to the terms and conditions of the BEI Resources Material Transfer Agreement (MTA). The MTA is available on our Web site at [www.beiresources.org](http://www.beiresources.org).

While BEI Resources uses reasonable efforts to include accurate and up-to-date information on this product sheet, neither ATCC® nor the U.S. Government makes any warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. Neither ATCC® nor the U.S. Government warrants that such information has been confirmed to be accurate.

This product is sent with the condition that you are responsible for its safe storage, handling, use and disposal. ATCC® and the U.S. Government are not liable for any damages or injuries arising from receipt and/or use of this product. While reasonable effort is made to ensure authenticity and reliability of materials on deposit, the U.S. Government, ATCC®, their suppliers and contributors to BEI Resources are not liable for damages arising from the misidentification or misrepresentation of products.

**Use Restrictions:**

**This material is distributed for internal research, non-commercial purposes only.** This material, its product or its derivatives may not be distributed to third parties. Except as performed under a U.S. Government contract, individuals contemplating commercial use of the material, its products or its derivatives must contact the contributor to determine if a license is required. U.S. Government contractors may need a license before first commercial sale.

**References:**

1. Bloom, J. and A. Balasz, Personal Communication.
2. Crawford, K. H. D., et al. “Dynamics of Neutralizing Antibody Titers in the Months after SARS-CoV-2 Infection.” *J. Infect. Dis.* 223 (2021): 197-205. PubMed: 33000143.

3. Crawford, K. H. D., et al. "Protocol and Reagents for Pseudotyping Lentiviral Particles with SARS-CoV-2 Spike Protein for Neutralization Assays." *Viruses* 12 (2020): E513. PubMed: 32384820.
4. Hulswit, R. J. G., C. A. M. de Haan and B. -J. Bosch. "Coronavirus Spike Protein and Tropism Changes." *Adv. Virus Res.* 96 (2016): 29-57. PubMed: 27712627.

ATCC® is a trademark of the American Type Culture Collection.



**APPENDIX I: NEUTRALIZATION ASSAY**

Reagents/Equipment		
Reagent/Equipment	Source	Catalog Number
HEK-293T-hACE2 Cell Line	BEI Resources	NR-52511
Polyclonal Anti-SARS-Related Coronavirus 2 Spike Glycoprotein (IgG, Rabbit)	BEI Resources	NR-52947
Polybrene	Sigma	TR1003-G
Lumitrac Plate	Greiner	655075

Note: Recommended maximum antibody dilution is 5 µg/mL per well, with 3-fold dilutions. Black walled, clear bottom or white-walled 96 well plates can be used for luciferase assays.

1. Plate HEK293T-hACE2 cells at a density of 1.5 x 10<sup>4</sup> cells per well in a 96 well plate.
2. In a separate U-bottom 96 well plate prepare antibody dilutions in cell culture medium (DMEM supplemented with 10% FBS) in a final volume of 30 µL.
3. Prepare control wells
  - a. Negative control (Mock)- add 50 µL cell culture medium
  - b. Positive control (Spike-LP)- add 30 µL cell culture medium + 20 µL working solution (step 4)
4. Add DMEM supplemented with 10% FBS to NR-53819 vial for a final volume of 2 mL (working solution); Dispense 20 µL working solution per well with antibody. Incubate plate 1 hour at 37°C under slow agitation.
5. Gently remove medium from cells; Add 10 µL of 30 µg/mL polybrene per well.
6. Inoculate cells with 50 µL Spike-LP/antibody mix from U-bottom plate per well; Incubate cells at 37°C with 5% CO<sub>2</sub> for 1 to 3 hours. Add 90 µL cell culture medium and incubate for 48 hours at 37°C with 5% CO<sub>2</sub>.
7. Perform luciferase assay or analyze GFP by microscopy.